



# Kangaroo by Jenx

## Technical Specifications



|                                              | Measurement    | Units          | Kangaroo DNR-55   |
|----------------------------------------------|----------------|----------------|-------------------|
| Age Range                                    | approx.        | <b>yrs</b>     | <b>3 - 7</b>      |
| User Weight                                  | max            | <b>kg</b>      | <b>35</b>         |
| Chest Height (footplate to top of chest pad) | min - max      | <b>mm</b>      | <b>740 - 1000</b> |
| Chest Pad Width                              | min - max      | <b>mm</b>      | <b>205 - 300</b>  |
| Hip Pad Width                                | min - max      | <b>degrees</b> | <b>205 - 300</b>  |
| Tray                                         | width x length | <b>mm</b>      | <b>580 x 425</b>  |
| Base Footprint                               | width x length | <b>mm</b>      | <b>585 x 705</b>  |