



Corner Seat by Jenx

Technical Specifications



Measurement		Units	Size 1	Size 2	Size 3
Age Range (approx)	approx.	mths - yrs	6 mths - 2 yrs	1 yrs - 6 yrs	4 yrs - 10 yrs
User Weight	max	kg	25	35	50
Total Back Support Height	min - max	mm	315 - 425	440 - 620	500 - 680
Lower Back Support Height	max	mm	165	235	235
Hip Width	max	mm	330	390	495
Footprint	width x length	mm	500 x 500	540 x 565	600 x 655