



Junior⁺ by Jenx

Technical Specifications



IMPORTANT INFORMATION: Please be aware, some minimum and maximum measurements are achieved by the use of additional accessory components.		Measurement	Units	Junior⁺ 1 JUC21	Junior⁺ 2 JUC22	Junior⁺ 3 JUC23
Age Range	approx.		mths - yrs	18mths - 5yrs	3yrs - 8yrs	5yrs - 12yrs
User Weight	max		kg	50	50	50
Back Angle	recline - forward		degrees	15° - 15°	15° - 15°	15° - 15°
Back Pad Height from Seat to Top	min - max		mm	345 - 435	345 - 435	410 - 500
Chest Width	min - max		mm	185 - 280	185 - 280	210 - 300
Hip Width	min - max		mm	150 - 375	150 - 375	150 - 375
Seat Depth	min - max		mm	210 - 305	240 - 355	260 - 440
Lower Leg Length	min - max		mm	195 - 270	235 - 350	250 - 405
Footplate Angle	dorsi - plantar		degrees	5° - 35°	5° - 35°	5° - 35°
X and Y Base Tilt-In-Space Angle	supine - prone		degrees	30° - 10°	30° - 10°	30° - 10°
X Base Footprint	width x length		mm	625 x 625	625 x 625	625 x 625
Y Base Footprint	width x length		mm	660 x 700	660 x 700	660 x 700
X Base Black Seat Height From Floor	min - max		mm	345 - 515	345 - 515	345 - 515
X Base White Seat Height From Floor	min - max		mm	395 - 595	395 - 595	395 - 595
Y Base Seat Height From Floor	min - max		mm	335 - 755	335 - 755	335 - 755