

Seating

Standing

Sleeping

Therapy

Kangaroo by Jenx

A flexible and adjustable standing system, providing assisted upright positioning for older children between 3 and 7 years with moderate ability

INNOVATIVE FEATURE

Adjustable Width Chest Side Supports

These supports are used to give extra side support to keep your child in mid-line and symmetrical. They can be adjusted for width to cater for larger or smaller children, and can be angled to give a really close fit.

INNOVATIVE FEATURE

Multi-Adjustable Knee Supports

The knee support can be adjusted for height, angle and position, forwards or backwards. The individual knee cups can be adjusted horizontally to achieve the ideal position for your child.



INNOVATIVE FEATURE

Adjustable Hip Support

The hip support can be adjusted for height and angle, forwards and backwards as well as width, perfect for adjusting to each particular child and also for children with feeding pegs.

INNOVATIVE FEATURE

Adjustable Height

The height adjustment on Kangaroo means that each child is fully supported as they grow throughout the recommended age range.



Visit jenx.com for more information



Kangaroo by Jenx

Technical Specifications

	Measurement	Units	Kangaroo DNR-55
Age Range	approx.	yrs	3 - 7
User Weight	max	kg	35
Chest Height (footplate to top of chest pad)	min - max	mm	740 - 1000
Chest Pad Width	min - max	mm	205 - 300
Hip Pad Width	min - max	degrees	205 - 300
Tray	width x length	mm	580 x 425
Base Footprint	width x length	mm	585 x 705

Find Your Distributor

To locate our carefully selected distributor for your region, please visit:

jenx.com/where-to-buy

