



Dreama[™] by Jenx

Technical Specifications



	Measurement	Units	Cot DC06	Cot Bed DC51	Single Bed DC03	Extra Long Single Bed DC15
Age Range	approx.	mths-yrs	9mths - 4yrs	9mths - 7yrs	9mths - Adult	9mths - Adult
User Weight	max	kg	80	100	110	110
Single Cushion Size	length x width x height	mm	560 x 95 x 150	700 x 95 x 150	900 x 95 x 150	900 x 95 x 150
Double Cushion Size	length x width x height	mm	560 x 195 x 150	700 x 195 x 150	900 x 195 x 150	900 x 195 x 150
No. of Cushions Supplied	(single:double)	mm	8:2	10:2	15:2	16:2
Height of System (Inc. Base)	max	mm	170	170	170	170
Footprint	width x length	mm	560 x 1200	700 x 1400	900 x 1900	900 x 2000